

Journal of Veterinary Diagnostic Investigation

<http://vdi.sagepub.com/>

Evaluation of the Mythic 18 hematology analyzer for use with canine, feline, and equine samples

Andrea K. Wassmuth, Barbara Riond, Regina Hofmann-Lehmann and Hans Lutz

J VET Diagn Invest 2011 23: 436

DOI: 10.1177/1040638711403416

The online version of this article can be found at:

<http://vdi.sagepub.com/content/23/3/436>

Published by:



<http://www.sagepublications.com>

On behalf of:



Official Publication of the American Association of Veterinary Laboratory Diagnosticians, Inc.

Additional services and information for *Journal of Veterinary Diagnostic Investigation* can be found at:

Email Alerts: <http://vdi.sagepub.com/cgi/alerts>

Subscriptions: <http://vdi.sagepub.com/subscriptions>

Reprints: <http://www.sagepub.com/journalsReprints.nav>

Permissions: <http://www.sagepub.com/journalsPermissions.nav>

>> [Version of Record](#) - May 1, 2011

[What is This?](#)

Evaluation of the Mythic 18 hematology analyzer for use with canine, feline, and equine samples

Andrea K. Wassmuth, Barbara Riond,¹ Regina Hofmann-Lehmann, Hans Lutz

Abstract. The Mythic 18 is a fully automated hematology bench-top analyzer using impedance technology for a complete blood cell count (CBC) and a 3-part white blood cell count (WBC) differential. The purpose of the current study was to evaluate the Mythic for assessment of agreement, precision, linearity, carry-over, stability, and usability under practice conditions. Ethylenediamine tetra-acetic acid–blood samples from 122 dogs, 140 cats, and 123 horses were analyzed with the Mythic and reference methods (Sysmex XT-2000iV, manual hematocrit, and microscopic WBC differentiation). Pearson’s coefficient of correlation, Passing–Bablok regression analysis, and Bland–Altman difference plots were performed to determine agreement. For precision, standard deviation and coefficients of variation were calculated. Linearity was determined according to Emancipator–Kroll. Red blood cell parameters showed excellent correlation and small biases, except for red cell distribution width and mean corpuscular hemoglobin concentration. Total WBC correlated excellently in canine and equine and very well in feline samples. In 23 feline specimens with platelet aggregates, the Mythic overestimated WBC. In all 3 species, absolute granulocyte counts correlated excellently. Equine lymphocyte counts showed good correlation whereas canine and feline lymphocyte counts correlated poorly. Feline platelets showed good correlation with a negative bias. The instrument showed good to excellent precision. The whole 3-part differential was found to be accurate in horses. In dogs and cats, absolute granulocyte counts were reliable. As with all impedance-based hematological instruments, evaluation of a blood smear is absolutely indicated to check for the presence of platelet aggregates, to verify WBC differentiation, and to identify possible abnormalities.

Key words: Cats; dogs; horses; impedance hematology analyzer; Mythic 18; 3-part differential.

Introduction

Hematological results provide important information on the patient’s state of health, disease history, and response to treatment.²³ The invention of the Coulter cell counter and cell volume analyzer in 1956 highly reduced time-consuming manual work by automating the counting and sizing of cells.⁹ Since then, several affordable, automated bench-top hematology analyzers have been developed for in-clinic use.³ Most of these analyzers are primarily designed for human blood. When analyzing nonhuman hematology specimens, it is essential that the selected instrument is designed and validated for multispecies analysis.²⁰

The Mythic 18^a is an impedance-based hematology instrument originally designed for human application. To make the instrument suitable for veterinary application, settings for feline, canine, and equine blood samples have been developed in the Clinical laboratory of the Vetsuisse-Faculty of the University of Zurich (Zurich, Switzerland). The evaluation was conducted to assure the quality of the newly designed animal settings. Information on imprecision and inaccuracy of a hematological instrument are extremely

valuable for the users. Each type of instrument should therefore be validated for each species before using results for clinical purpose.

The objective of the present study was to validate the Mythic 18 for use with blood samples from healthy and diseased cats, dogs, and horses. To this end agreement, precision, linearity, carry-over, and sample stability were determined. Biases were judged with respect to their clinical relevance.

Material and methods

Blood samples

Fresh ethylenediamine tetra-acetic acid (EDTA)-K₃ blood samples from 122 dogs, 140 cats, and 123 horses from the

From the Clinical Laboratory, Vetsuisse Faculty, University of Zurich, Switzerland.

¹Corresponding Author: Barbara Riond, Clinical Laboratory, Vetsuisse-Faculty, University of Zurich, Winterthurerstrasse 260, 8057 Zurich, Switzerland. briond@vetclinics.uzh.ch

Small Animal Clinic and the Clinic for Horses, at the Vetsuisse-Faculty, University of Zurich were analyzed on the Sysmex XT-2000iV^b and reference methods, and usually with a time delay of 1.5 hr (until the routine work was finished) on the Mythic 18. All samples were collected by venipuncture regardless of sex, age, or breed and sent to the clinical laboratory in the framework of routine work to check the health status. Sample collection took place between May and December 2009. Complete sample analysis was performed within 6 hr after collection, most of them within 4 hr. The aforementioned blood samples were used to assess agreement and precision. To determine the range of linear measurement, 2 blood samples from cats, 2 from dogs, and 1 from a horse were used. Additionally, platelet-enriched plasma from a horse was used to assess linearity of the platelet count. Carry-over of blood from one sample to the following sample, thereby checking the effectiveness of cleaning of the instrument, was assessed for each species using 2 EDTA blood samples. To determine the effect of aging of samples, blood samples from 6 dogs and cats and 8 horses were used.

Instruments and methods used

Mythic 18. The Mythic 18 is a fully automated in-house hematology analyzer performing hematological analyses on EDTA anticoagulated blood. Nineteen species profiles can be created. For counting the cellular blood components, the Mythic 18 uses the impedance technique. A cyanide-free spectrophotometry method is used to measure hemoglobin by formation of oxyhemoglobin at 555 nm. Hematocrit is measured by volume integration. The sample volume is 10 μ l. The instrument can determine 16 parameters in the normal mode and 18 in the research mode: white blood cell count (WBC) with absolute number and percentage of lymphocytes (LYM), monocytes (MONO), and granulocytes (GRAN), red blood cell count (RBC), hemoglobin concentration (HGB), hematocrit value (HCT), mean corpuscular volume (MCV), mean corpuscular hemoglobin (MCH), mean corpuscular hemoglobin concentration (MCHC), red cell distribution width (RDW), platelets (PLT), mean platelet volume (MPV), and for research, platelet crit (PCT) and platelet distribution width (PDW). For platelet counting, a floating threshold is used, whereas for RBC and WBC counts, the thresholds are predefined. Results are provided within 1 min on the liquid crystal display, printed out on a printer, and stored in the resident memory or in a USB key. Results were presented with flags; optionally, reference ranges can be reported. Additionally, the Mythic 18 shows histograms for WBC, RBC, and PLT. Prior to analysis, patient's data can be entered manually or with a barcode reader. The instrument also displays message codes and histogram flags. However, message codes and flags have not been adapted yet to feline, canine, and equine blood samples. Therefore, conclusions about the usefulness of these message codes and flags cannot be drawn at this time.

Mythic 18 provides a 3-part WBC differential in samples with WBC within a range of $0.9 \times 10^3/\mu$ l to $150 \times 10^3/\mu$ l. Quality control samples^c are supplied as blood samples with 3 levels of RBC, WBC, and PLT levels. Results of each lot can be viewed on the display of the instrument in tables and Levey–Jennings graphs. The instrument uses 3 reagents^d: a diluent, a lysis reagent, and a cleaning solution.

Sysmex XT-2000iV. The Sysmex XT-2000iV, equipped with software version 10b, was used as the reference instrument for total WBC, WBC differentiation, RBC, RBC indices, HGB, RDW, PLT count, and MPV. It is a fully automated hematology analyzer for animal blood providing 30 parameters. The impedance method with hydrodynamic focusing is used for RBC (RBC-I), HCT, and PLT (PLT-I). With these results, MCV, MCH, MCHC, RDW, MPV, and PDW are calculated. A flow cytometry device based on a sheath-flow and a semiconductor laser is used as an optical method for platelets (PLT-O) in cats, WBC, and WBC differentiation. Hemoglobin concentration is measured spectrophotometrically with a cyanide-free (sodium lauryl sulfate) method.

Manual methods. Manual HCT measurement was done with microhematocrit capillary tubes centrifuged at $13,000 \times g$ for 5 min in a microhematocrit centrifuge.¹⁰ Blood smears were stained using an automated staining instrument.^c Microscopic differentiation of 2 modified Wright–stained blood smears, 100 WBC each, was conducted by 2 technicians with 10 years each of veterinary hematology experience. The results were used to calculate the absolute number of LYM, MONO, and GRAN counts by multiplying the percentage from the 200-cell count of each cell type with the total WBC from the Sysmex XT-2000iV.

Assessment of agreement

Analytical accuracy is defined by the International Council for Standardization in Hematology as a measure of agreement between the measured value of an analyte and its “true” value.⁷ To determine accuracy, agreement between the results of the evaluated instrument and the results of a reference instrument were investigated. In the present study, assessment of agreement was determined by comparing the results of the Mythic 18 with those of the reference instrument, the manual HCT, and the microscopic differentiation. Sysmex XT-2000iV is widely used and accepted in veterinary clinical laboratories, and validation studies were conducted on the Sysmex XT-2000iV for cats²¹ and cats, dogs, and horses.^{13,14} For comparing results of granulocytes of the Mythic 18, results of the neutrophils, eosinophils, and basophils of the reference methods were added.

Precision and linearity

Within-series precision of the instrument was determined for each of the investigated species for low, normal, and high

WBC values based on multiple analyses (more than 12 consecutive times). During the analysis, the sample was gently mixed. Afterwards, mean, standard deviation, and coefficient of variation as a measurement of the random error were calculated for all parameters. Precision from day to day was measured using commercially available quality control blood of low, intermediate, and high levels, which were analyzed once daily prior to analyzing patients samples over a 20-day period.

The linearity of the measurement range was assessed in all 4 species to determine the analytical range. Mythic 18 has a reportable range for WBC ($-150 \times 10^3/\mu\text{l}$), RBC ($-15 \times 10^6/\mu\text{l}$), HCT (-72%), and PLT ($-4.000 \times 10^3/\mu\text{l}$). The linearity of the measurement range was determined for WBC, RBC, HCT, HGB, and PLT by analyzing a series dilution of EDTA-K₃ anticoagulated blood in triplicate. For cat and dog, 2 blood samples were used, 1 with high WBC to determine WBC linearity (5 ml) and 1 (cat, 12 ml; dog, 10 ml) for the remaining parameters. One equine sample (20 ml of EDTA-blood) was used for RBC, HGB, and HCT; additionally, platelet-enriched plasma of a horse was used to determine PLT linearity. The blood samples were centrifuged^f at $390 \times g$ for 10 min to receive results below and above the reference range. The plasma was then removed from the blood cells. Afterwards, concentrated blood cells were diluted with 0.9% saline solution in steps of 10%, to achieve a dilution series of 0–100% blood cell concentrate.

Carry-over

Carry-over was studied to assess if transfer of blood from a single sample will cause a falsely higher result in the following sample. For each species, 2 patient samples, 1 with high WBC, were analyzed twice followed by 3 replicates of diluents.¹²

Cell aging

Cell aging studies were performed with blood samples from 6 cats, 6 dogs, and 8 samples from horses. They were analyzed at time point 1, 2, 4, 6, 8, 24, 32, and 48 hr after collection to calculate stability. The blood was stored at room temperature during the whole experiment. For each parameter, the difference in mean of results between each analysis and time point 1 hr was calculated. In cats, only RBC parameters were investigated.

Statistical analysis

All data were entered manually in a Microsoft Excel^g spreadsheet. The Microsoft Excel add-in Analyse-it^h was used for statistical analyses. For each parameter and each investigated species, Pearson coefficient of correlation (r), linear regression analysis according to Passing–Bablok provided intercept and slope with 95% confidence interval,² and

Table 1. Reference values of hematological parameters for cats, dogs, and horses used in the current study (lower 5% reference limit and upper 95% reference limit).

Parameter*	Cat	Dog	Horse
WBC ($\times 10^3/\mu\text{l}$)	4.6–12.8	4.7–11.3	4.7–8.2
LYM ($/\mu\text{l}$)	1,050–6,000	1,154–3,399	1,020–3,472
MONO ($/\mu\text{l}$)	46–678	198–917	0–184
Bands ($/\mu\text{l}$)	0–123	0–84	0–75
Neutrophils ($/\mu\text{l}$)	2,315–10,011	2,496–7,437	3,021–5,775
Eosinophils ($/\mu\text{l}$)	100–600	119–1,287	0–216
Basophils ($/\mu\text{l}$)	0–143	0–82	0–66
RBC ($\times 10^6/\mu\text{l}$)	7–10.7	6.1–8.1	6.2–9
HGB (g/dl)	11.3–15.5	14.4–19.1	10.8–14.9
HCT (%)	33–45	42–55	30–42
MCH (pg)	14–17	23–26	15–18
MCV (fl)	41–49	64–73	41–50
MCHC (g/dl)	33–36	34–36	35–37
PLT ($\times 10^3/\mu\text{l}$)	180–680	130–394	119–250

*WBC = white blood cell count; LYM = absolute number of lymphocytes; MONO = monocytes; RBC = red blood cell count; HGB = hemoglobin concentration; HCT = hematocrit value; MCH = mean corpuscular hemoglobin; MCV = mean corpuscular volume; MCHC = mean corpuscular hemoglobin concentration; PLT = platelets.

Bland–Altman difference plot with biases and 95% limits of agreement¹ were calculated. Coefficient of correlation was considered excellent if $r \geq 0.95$, very good if $r = 0.90$ – 0.94 , good if $r = 0.80$ – 0.89 , fair if $r = 0.59$ – 0.79 , and poor if $r < 0.59$.²² For precision analysis, standard deviation (SD) and coefficient of variation (CV) were calculated for each level, parameter, and species. The degree of linearity was determined with Analyse-it according to Emancipator–Kroll.¹¹

For WBC, RBC, HGB, and PLT, the percentage of carry-over was calculated using the following formula:

$$\% \text{ carry-over} = \left\{ \frac{[(\text{result empty cycle 1}) - (\text{result empty cycle 3})]}{[(\text{sample 2}) - (\text{result empty cycle 3})]} \right\} \times 100.$$

The stability of the blood samples was reviewed for statistically significant changes using the Friedman test and the Dunn multiple comparison post-test with GraphPad Prism.¹ Statistical significance was tested for the result of the first hour compared with the results of the following time points. Statistical significance was defined as P value < 0.05 .

Clinical relevance

For each sample, the data from the Mythic 18 and the reference methods were compared with hematology reference values that were established and currently used in the Clinical Laboratory of the Vetsuisse-Faculty University of Zurich (Table 1). The results were judged to be below or above the

Table 2. Agreement between results of the Mythic 18 hematology analyzer^a and results of the reference methods (Sysmex XT-2000iV hematology analyzer^b and microhematocrit tubes).

Parameter*	Species	Coefficient of correlation	Intercept [†]	Slope [†]	Bias [‡]	No. of samples
WBC	Cat	0.94	0.26 (-0.14–0.61)	0.91 (0.88–0.95)	-0.072 (-6.959–6.815)	129
	Dog	0.99	0.98 (0.68–1.36)	0.95 (0.92–0.98)	0.229 (-3.651–4.110)	122
	Horse	0.98	0.38 (0.17–0.56)	0.94 (0.92–0.97)	-0.126 (-1.496–1.243)	123
RBC	Cat	0.99	0.51 (0.37–0.66)	0.95 (0.93–0.97)	0.090 (-0.400–0.581)	138
	Dog	0.99	0.26 (0.15–0.40)	1.00 (0.98–1.02)	0.241 (-0.096–0.578)	122
	Horse	0.98	0.38 (0.21–0.57)	0.93 (0.91–0.96)	-0.140 (-0.693–0.414)	123
HGB	Cat	0.99	0.38 (0.21–0.54)	0.92 (0.90–0.93)	-0.511 (-1.256–0.234)	138
	Dog	1.00	1.02 (0.81–1.26)	0.93 (0.92–0.95)	0.10 (-0.57–0.76)	122
	Horse	0.98	0.48 (0.08–0.76)	0.94 (0.92–0.98)	-0.25 (-1.09–0.58)	123
HCT	Cat	0.99	2.05 (1.23–2.74)	0.94 (0.92–0.97)	0.16 (-2.15–2.48)	135
	Dog	0.99	1.20 (0.05–2.13)	0.95 (0.93–0.98)	-0.79 (-3.42–1.83)	121
	Horse	0.99	1.36 (0.43–2.13)	0.96 (0.93–0.98)	-0.17 (-1.95–1.61)	123
MCV	Cat	0.95	4.65 (2.50–7.12)	0.91 (0.85–0.96)	0.86 (-2.62–4.33)	138
	Dog	0.96	8.01 (4.54–11.37)	0.83 (0.78–0.88)	-3.16 (-6.03–0.28)	122
	Horse	0.94	3.62 (1.08–5.82)	0.90 (0.85–0.95)	-1.02 (-3.68–1.63)	123
MCH	Cat	0.97	-0.08 (-0.80–0.49)	0.95 (0.91–1.00)	-0.84 (-1.62–0.06)	138
	Dog	0.92	-3.01 (-5.25–0.80)	1.10 (1.00–1.20)	-0.65 (-1.91–0.62)	122
	Horse	0.95	0.00 (-1.47–0.00)	1.00 (1.00–1.09)	-0.02 (-0.72–0.69)	123
MCHC	Cat	0.67	8.58 (4.03–12.54)	0.67 (0.55–0.80)	-2.79 (-5.80–0.21)	138
	Dog	0.65	0.50 (-8.63–7.24)	1.00 (0.80–1.26)	0.65 (-1.80–3.11)	122
	Horse	0.50	8.42 (0.90–14.53)	0.79 (0.62–1.00)	0.75 (-1.51–3.02)	123
RDW	Cat	0.37	7.55 (4.19–10.17)	0.45 (0.33–0.59)	-5.19 (-10.24–0.15)	138
	Dog	0.73	4.54 (2.78–5.89)	0.56 (0.47–0.68)	-2.33 (-4.92–0.27)	121
	Horse	0.47	12.39 (10.00–14.41)	0.28 (0.19–0.38)	-4.98 (-8.83–1.13)	123
PLT	Cat	0.80	-9.47 (-52.65–23.06)	0.88 (0.74–1.05)	-38.3 (-225.5–149.0)	90
	Dog	0.97	-8.06 (-27.00–7.22)	1.15 (1.10–1.22)	42.5 (-73.9–158.8)	121
	Horse	0.84	-18.08 (-37.57–4.82)	1.04 (0.93–1.16)	1.3 (-82.3–84.9)	117
MPV	Dog	0.73	0.72 (-0.08–1.44)	0.69 (0.62–0.77)	-2.25 (-3.77–0.73)	106
	Horse	0.80	0.10 (-1.37–0.10)	1.00 (1.00–1.20)	0.10 (-0.50–0.69)	100

*WBC = white blood cell count; RBC = red blood cell count; HGB = hemoglobin concentration; HCT = hematocrit value; MCV = mean corpuscular volume; MCH = mean corpuscular hemoglobin; MCHC = mean corpuscular hemoglobin concentration; RDW = red cell distribution width; MPV = mean platelet volume.

[†]Numbers in parentheses are 95% confidence intervals.

[‡]Numbers in parentheses are 95% limits of agreement.

reference range, and the resulting interpretations from the Mythic 18 and the reference methods were compiled and compared to each other.

Results

Assessment of agreement

Pearson coefficient of correlation, intercept, and slope with 95% confidence intervals (CI) calculated by Passing–Bablok regression analysis, and biases with their 95% limits of agreement calculated by Bland–Altman difference plot are presented in Tables 2 and 3. Table 2 shows results for WBC, RBC, and PLT, and Table 3 presents results of the WBC differentiation compared with results from the Sysmex XT-2000iV and results of the manual WBC differentiation. Furthermore, linear regression analysis by Passing–Bablok

and Bland–Altman difference plots are presented (Fig. 1 [cat], Fig. 2 [dog], and Fig. 3 [horse]). The results of the Mythic 18 for RBC, HGB concentration, HCT, and WBC (except for the cat) showed excellent correlation with the results provided by the reference instrument Sysmex XT-2000iV and manual HCT. The WBC results from the Mythic 18 were compared with the optical WBC results of the Sysmex XT-2000iV. Systematic errors with very small biases were observed in all 3 species for WBC. The RBC showed an excellent result for dogs with a small bias due to a constant systemic error. For HGB levels, a small proportional systemic error was seen in all investigated species. The MCV values showed excellent correlation in cat and dog and very good correlation in horses, with a systemic error and negative biases for horses and dogs. The feline MCH showed an excellent correlation with a small negative bias due to a

Table 3. Agreement between the WBC differentiation for absolute numbers of the Mythic 18 hematology analyzer^a compared with results from the Sysmex XT-2000iV hematology analyzer^b and manual differentiation.

Parameter	Species	Coefficient of correlation	Intercept*	Slope*	Bias [†]	No. of samples	
Lymphocyte	Sysmex	Cat	0.52	0.09 (-0.51-0.54)	1.27 (1.00-1.63)	1.014 (-3.324-5.351)	115
		Dog	0.19	0.25 (-0.56-0.69)	1.11 (0.82-1.62)	0.687 (-2.683-4.056)	118
		Horse	0.89	0.16 (0.00-0.32)	0.96 (0.89-1.06)	0.127 (-0.957-1.212)	122
	Manual	Cat	0.49	0.12 (-0.28-0.59)	1.50 (1.15-1.90)	1.348 (-3.280-5.975)	120
		Dog	0.14	0.42 (-0.25-0.91)	1.24 (0.88-1.86)	1.095 (-2.435-4.625)	117
		Horse	0.87	0.37 (0.16-0.55)	0.88 (0.80-0.97)	0.139 (-1.129-1.408)	123
Monocyte	Sysmex	Cat	0.37	-0.05 (-0.19-0.05)	1.74 (1.33-2.31)	0.198 (-0.902-1.299)	117
		Dog	0.63	0.29 (0.19-0.38)	0.70 (0.57-0.86)	0.010 (-1.158-1.178)	119
		Horse	0.45	0.02 (-0.06-0.09)	0.77 (0.56-1.00)	-0.065 (-0.468-0.339)	122
	Manual	Cat	0.60	0.04 (-0.06-0.13)	1.43 (1.11-1.82)	0.253 (-0.860-1.366)	120
		Dog	0.57	0.34 (0.24-0.45)	0.63 (0.45-0.81)	0.047 (-1.366-1.460)	117
		Horse	0.24	0.15 (0.10-0.18)	0.61 (0.42-0.83)	0.074 (-0.402-0.550)	123
Granulocyte	Sysmex	Cat	0.97	0.52 (0.25-0.73)	0.83 (0.79-0.87)	-1.352 (-7.149-4.445)	115
		Dog	0.99	1.01 (0.72-1.29)	0.89 (0.86-0.93)	-0.451 (-5.235-4.332)	117
		Horse	0.98	0.24 (0.04-0.45)	0.93 (0.89-0.97)	-0.203 (-1.641-1.235)	122
	Manual	Cat	0.97	0.47 (0.19-0.82)	0.81 (0.77-0.85)	-1.643 (-7.967-4.682)	120
		Dog	0.99	1.05 (0.72-1.31)	0.86 (0.83-0.90)	-0.997 (-6.441-4.446)	117
		Horse	0.97	0.23 (0.03-0.52)	0.91 (0.86-0.95)	-0.347 (-2.013-1.318)	123

*Numbers in parentheses are 95% confidence intervals.

[†]Numbers in parentheses are 95% limits of agreement.

constant systematic error, whereas the dog showed a proportional systemic error with a negative bias. For MCHC and PLT counts, a proportional systemic error was seen in all 3 species. In dogs, the Mythic 18 overestimated high PLT counts compared to the reference instrument. The MPV results for feline samples were not available, because the Sysmex XT-2000iV determines PLT counts optically via flow cytometry. The 3-part WBC differential showed for GRAN counts (absolute numbers) the best correlation and the smallest bias in all 3 species. The LYM counts showed a strong positive bias in cats and dogs with wide 95% limits of agreement. In horses, correlation was found to be good with a small bias. Results for MONO counts showed only fair correlation in canine samples and poor correlation in feline and equine samples. Some results of the WBC differential of the Sysmex XT-2000iV were excluded from statistical analysis due to the inability of the Sysmex XT-2000iV to differentiate WBC: 1 equine, 2 canine, and 4 feline blood samples. In the equine sample, both canine samples, and 3 of the 4 feline samples, the Sysmex XT-2000iV misclassified a left shift. In the remaining feline sample, normoblasts (52 normoblasts per 100 WBC) were seen in the blood smear while the Sysmex XT-2000iV classified them falsely to LYM.

Precision and linearity


Mean values and CVs from the precision study are presented in Table 4 (precision in series) and Table 5 (precision from day to day). Results of the linearity study are presented in Table 6. Linearity plots are shown in Figures 4 (cat), 5 (dog), and 6 (horse). For all tested parameters, the instrument demonstrated good linearity. The tested ranges of linearity were within the ranges provided by the manufacturer for human blood except for canine WBC.

Carry-over

Table 7 presents the results of the carry-over experiment ("value" in Table 7 represents the result of the second sample analysis). The results of the second and third diluent analyses were always 0. All results for carry-over lie in the range provided by the manufacturer (<1%), except 1 sample for feline WBC.

Cell aging

Table 8 shows the results of the cell aging study. The first significant changes appeared after 6 hr for HGB in the feline







Figure 1. Bland–Altman analyses respectively Passing–Bablok regression for feline accuracy results. Comparison of the Mythic 18 hematology analyzer^a with the Sysmex XT-2000iV hematology analyzer^b resp. manual hematocrit. For feline white blood cell count (WBC; **A**); absolute number of lymphocytes (LYM #; **B**), monocytes (MONO #; **C**), granulocytes (GRAN #; **D**); red blood cell count (RBC; **E**); hemoglobin concentration (HGB; **F**); hematocrit value (HCT; **G**); mean corpuscular volume (MCV; **H**); mean corpuscular hemoglobin (MCH; **I**); mean corpuscular hemoglobin concentration (MCHC; **J**); red cell distribution width (RDW; **K**); and platelets (PLT; **L**); In the Passing–Bablok regression plots, the thin gray line is the line of identity ($y = x$), and the thick black is the line of best fit. In Bland–Altman difference plots, the thin horizontal line (0 at the y-axis) is the line of identity, and the thick black line indicates the bias (mean difference between methods), with their confidence intervals as thin dashed lines. The thick dashed horizontal lines are the 95% limits of agreement with their 95% confidence intervals.



(continued)





Figure 2. Bland–Altman analyses respectively Passing–Bablok regression for canine accuracy results. Comparison of the Mythic 18 hematology analyzer^a with the Sysmex XT-2000iV hematology analyzer^b resp. manual hematocrit. For canine white blood cell count (WBC; **A**); absolute number of lymphocytes (LYM #; **B**), monocytes (MONO #; **C**), granulocytes (GRAN #; **D**); red blood cell count (RBC; **E**); hemoglobin concentration (HGB; **F**); hematocrit value (HCT; **G**); mean corpuscular volume (MCV; **H**); mean corpuscular hemoglobin (MCH; **I**); mean corpuscular hemoglobin concentration (MCHC; **J**); red cell distribution width (RDW; **K**); platelets (PLT; **L**); and mean platelet volume (MPV; **M**). In the Passing–Bablok regression plots, the thin gray line is the line of identity ($y = x$), and the thick black line is the line of best fit. In Bland–Altman difference plots, the thin horizontal line (0 at the y-axis) is the line of identity, and the thick black line indicates the bias (mean difference between methods), with their confidence intervals as thin dashed lines. The thick dashed horizontal lines are the 95% limits of agreement with their 95% confidence intervals.



(continued)




Figure 3. Bland–Altman analyses respectively Passing–Bablok regression for equine accuracy results. Comparison of the Mythic 18 hematology analyzer^a with the Sysmex XT-2000iV hematology analyzer^b resp. manual hematocrit. For equine white blood cell count (WBC; **A**); absolute number of lymphocytes (LYM #; **B**), monocytes (MONO #; **C**), granulocytes (GRAN #; **D**); red blood cell count (RBC; **E**); hemoglobin concentration (HGB; **F**); hematocrit value (HCT; **G**); mean corpuscular volume (MCV; **H**); mean corpuscular hemoglobin (MCH; **I**); mean corpuscular hemoglobin concentration (MCHC; **J**); red cell distribution width (RDW; **K**); platelets (PLT; **L**); and mean platelet volume (MPV; **M**); In the Passing–Bablok regression plots, the thin gray line is the line of identity ($y = x$), and the thick black line is the line of best fit. In Bland–Altman difference plots, the thin horizontal line (0 at the y-axis) is the line of identity, and the thick black line indicates the bias (mean difference between methods), with their confidence intervals as thin dashed lines. The thick dashed horizontal lines are the 95% limits of agreement with their 95% confidence intervals.

samples, but not at the remaining time points. In canine blood samples, a significant change in RBC was detected after 24 hr, but not at the remaining time points. Feline MCHC showed significant changes (i.e., a decrease) at time points 10, 24, 32, and 48 hr; in horses, such changes were also found at time points 32 hr and 48 hr. After 24 hr, equine WBC values started to decrease significantly, and the LYM-GRAN ratio moved in favor of LYM count. Canine and feline samples presented significant changes (i.e., an increase) for MCV values after 24, 32, and 48 hr. At time point 48 hr, feline HCT values, canine MCH values, and canine MONO% started to increase statistically significant.

Clinical relevance

Some of the results deviate from those determined by the reference methods. In Table 9, the number of results that deviate and their clinical relevance are compiled. Table 10 presents the abnormalities seen in the blood smear during manual WBC differentiation.

Discussion

The current study is based on the comparison of the 3-part WBC differential of the Mythic 18 with both manual and electronic WBC differentiation. Manual chamber counting of WBC, RBC, and PLT are well known as gold standard techniques. However, these methods show high imprecision due to the limited quantity of counted cells, artifacts, and classification of the cells.^{8,10,14} Therefore, electronic blood cell counting methods have mainly replaced the former gold standard techniques. Microscopic WBC differentiation of a blood smear is still mandatory to confirm WBC abnormalities and to rule out the presence of platelet clumps, RBC parasites, and blood cell precursors. Nevertheless, this manual technique is prone to high imprecision, especially for those cells that are underrepresented in the blood.²¹

Results for RBC parameter of the Mythic 18 showed very good to excellent agreement with the Sysmex XT-2000iV results, except for MCHC and RDW. Canine MCV values of the Mythic 18 showed a small negative bias compared to the

Table 4. Precision in series: mean values and coefficients of variation for blood samples from cats, dogs, and horses with low, normal, and high values for total white blood cell count.*

Parameter	Species	Low		Normal		High	
		Mean	CV	Mean	CV	Mean	CV
WBC ($\times 10^3/\mu\text{l}$)	Cat	1.96	2.50	6.81	1.50	40.92	1.17
	Dog	2.99	1.92	6.77	1.68	80.71	0.98
	Horse	1.72	3.24	7.85	1.58	19.76	0.95
LYM ($\times 10^3/\mu\text{l}$)	Cat	1.21	4.75	1.87	3.06	7.08	6.12
	Dog	1.05	4.74	1.64	5.81	9.75	3.93
	Horse	0.69	9.32	1.96	4.13	3.27	4.07
MONO ($\times 10^3/\mu\text{l}$)	Cat	0.13	34.91	0.2	0	3.58	10.16
	Dog	0.2	0	0.8	6.45	3.12	4.25
	Horse	0.09	27.74	0.21	12.43	0.46	13.14
GRAN ($\times 10^3/\mu\text{l}$)	Cat	0.6	8.61	4.75	1.7	30.25	3.93
	Dog	1.77	3.23	4.33	4.17	67.84	1.34
	Horse	0.99	2.59	5.66	1.84	16.05	1.88
RBC ($\times 10^6/\mu\text{l}$)	Cat	6.96	1.30	9.17	0.72	6.04	0.94
	Dog	3.99	0.75	7.85	0.54	4.87	1.32
	Horse	6.25	1.18	8.85	1.01	11.48	0.78
HGB (g/dl)	Cat	8.28	0.79	12.91	0.53	7.09	0.62
	Dog	8.65	0.93	15.79	0.56	9.99	1.03
	Horse	11.49	1.13	15.06	0.85	19.73	0.45
HCT (%)	Cat	26.38	1.66	39.51	0.82	26.21	1.02
	Dog	26.50	0.86	46.36	0.68	31.14	1.40
	Horse	30.17	1.47	40.43	0.95	51.91	1.13
MCV (fl)	Cat	37.93	0.67	43.09	0.31	43.39	0.58
	Dog	66.37	0.34	59.05	0.43	63.92	0.43
	Horse	48.25	0.70	45.66	0.28	45.20	0.49
MCH (pg)	Cat	11.90	1.15	14.07	1.02	11.74	0.81
	Dog	21.68	0.97	20.11	0.83	20.49	1.08
	Horse	18.36	0.86	17.02	0.87	17.21	0.60
MCHC (g/dl)	Cat	31.40	1.36	32.67	1.07	27.06	1.05
	Dog	32.67	1.19	34.06	1.06	32.08	1.32
	Horse	38.09	1.37	37.26	0.90	38.01	0.96
RDW (%)	Cat	21.93	1.74	16.88	1.90	16.99	2.10
	Dog	16.32	2.60	16.06	2.27	13.23	1.60
	Horse	18.61	2.13	18.46	2.81	18.68	1.52
PLT ($\times 10^3/\mu\text{l}$)	Cat	133.47	5.65	249.73	4.63	170.07	6.91
	Dog	305.73	3.17	215.47	3.98	144.53	6.23
	Horse	105.79	10.24	205.36	4.67	73.21	7.66
MPV (fl)	Cat	9.01	2.68	9.65	0.99	9.72	2.10
	Dog	9.38	1.24	7.50	1.89	10.12	2.11
	Horse	8.01	1.97	7.69	2.08	6.94	3.52

*CV = coefficient of variation; WBC = white blood cell count; LYM = lymphocytes; MONO = monocytes; GRAN = granulocytes; RBC = red blood cell count; HGB = hemoglobin concentration; HCT = hematocrit value; MCV = mean corpuscular volume; MCH = mean corpuscular hemoglobin; MCHC = mean corpuscular hemoglobin concentration; RDW = red cell distribution width; PLT = platelets; MPV = mean platelet volume.

Sysmex XT-2000iV, which could lead to different clinical conclusion in some cases. Small changes of the HCT correction factor in the canine settings could improve MCV agreement between the Sysmex XT-2000iV and the Mythic 18. Otherwise, adjustment of the reference limits for canine

MCV would be indicated. Difference in the osmolarity of the diluents between the Sysmex XT-2000iV (250 mosm/kg) and the Mythic 18 (332 mosm/kg) can possibly cause a negative bias in canine MCV values. The hypotonic diluent of the Sysmex XT-2000iV causes swelling of the RBC whereas the

Table 5. Precision from day to day: mean value and coefficient of variation for control blood samples.*

Parameter	Low		Middle		High	
	Mean	CV	Mean	CV	Mean	CV
WBC ($\times 10^3/\mu\text{l}$)	2	3.2	7.4	2	19	0.8
LYM ($\times 10^3/\mu\text{l}$)	1.1	5.7	2.1	5.1	3.1	6.6
LYM (%)	55.1	3.3	29.1	3.7	16	5.8
MONO ($\times 10^3/\mu\text{l}$)	0.2	20.1	0.4	9.4	0.6	3.6
MONO (%)	11.7	8.5	5.8	5.7	3.2	5
GRAN ($\times 10^3/\mu\text{l}$)	0.7	7.9	4.8	2	15.4	1.2
GRAN (%)	33.2	3.7	65.1	1.7	80.8	1.3
RBC ($\times 10^6/\mu\text{l}$)	2.57	1.9	4.97	1.6	5.99	1.2
HGB (g/dl)	6.6	2.2	14.4	1.5	18.8	1.4
HCT (%)	17	2	36.9	1.2	48	1
MCV (fl)	66.2	0.8	74.2	0.7	80.2	0.6
MCH (pg)	25.6	1.9	28.9	1.6	31.4	1.3
MCHC (g/dl)	38.7	1.9	39	1.4	39.1	1.3
RDW (%)	16.8	3.4	16.2	2.9	14.2	2.7
PLT ($\times 10^3/\mu\text{l}$)	84	8.4	236	5.9	482	4.6
MPV (fl)	8.4	3.7	8	2.3	7.8	2

*CV = coefficient of variation; WBC = white blood cell count; LYM = lymphocytes; MONO = monocytes; GRAN = granulocytes; RBC = red blood cell count; HGB = hemoglobin concentration; HCT = hematocrit value; MCV = mean corpuscular volume; MCH = mean corpuscular hemoglobin; MCHC = mean corpuscular hemoglobin concentration; RDW = red cell distribution width; PLT = platelets; MPV = mean platelet volume.

relatively isotonic diluent of the Mythic 18 produces comparatively lower MCV values.⁴ The agreement for MCHC values is less satisfactory in all 3 evaluated species. Low correlation for this parameter has been reported in previous studies^{18,21,23} and can be mainly explained by the narrow concentration range of this parameter.

Generally, the Mythic 18 underestimated high total WBC in all 3 species on average by a few percentage points. In the feline samples, total WBC correlation is very good. However, in 23 out of 129 feline samples, the WBC were on average more than 1,000 WBC/ μl higher than those of the reference method in which the WBC are determined by an optical flow cytometry principle⁹ (Fig. 1A). This can readily be explained by the fact that feline PLT have a high tendency to aggregate. When 2 samples with the extremely high overestimation of more than 15,000 WBC/ μl were removed from the statistic as outliers, the coefficient of correlation improved to 0.97 and the bias decreased from -0.043 to -0.364. In all 23 samples, platelet aggregates could be identified in the blood smear. It is a well-known phenomenon in cats that platelet clumps or large platelets can cause falsely increased WBC and decreased PLT counts in impedance-based hematological instruments.^{10,16}

In horses, cats, and dogs, GRAN are predominant in the blood (Table 1), therefore imprecision for this leukocyte

Table 6. Range of linearity for canine and feline white blood cell count, red blood cell count parameter, and equine platelets in cats, dogs, and horses and for human blood samples.

Species	Parameter*	Range of linearity
Cat	WBC	-100 ($\times 10^3/\mu\text{l}$)
	RBC	-14 ($\times 10^6/\mu\text{l}$)
	HGB	-23.2 (g/dl)
	HCT	-70 (%)
Dog	WBC	-95 ($\times 10^3/\mu\text{l}$)
	RBC	-11.4 ($\times 10^6/\mu\text{l}$)
	HGB	-24.3 (g/dl)
Horse	HCT	-71 (%)
	RBC	-12.6 ($\times 10^6/\mu\text{l}$)
	HGB	-24.4 (g/dl)
	HCT	-62 (%)
Human blood [‡]	PLT [†]	-1.080 ($\times 10^3/\mu\text{l}$)
	WBC	0-100 ($\times 10^3/\mu\text{l}$)
	RBC	0.1-8 ($\times 10^6/\mu\text{l}$)
	HGB	0.5-24 (g/dl)
	HCT	5-70 (%)
	PLT	5-2,000 ($\times 10^3/\mu\text{l}$)

*WBC = white blood cell count; RBC = red blood cell count; HGB = hemoglobin concentration; HCT = hematocrit value; PLT = platelets.

[†]For this study, platelet-enriched plasma from a horse was used.

[‡]Human blood samples were obtained from Orphée SA, Geneva, Switzerland.

subtype is low, and the Mythic 18 showed excellent agreement with both the Sysmex XT-2000iV and the manual WBC differentiation.¹⁷ Agreement in canine and feline LYM counts in the Mythic 18 was not satisfactory. This finding has already been demonstrated for both species in the VetScan HMT,⁶ and for canine samples in the CA530-Vet¹⁷ and the Heska CBC (Becker M: 2007, A comparative study of seven in-house and two laboratory hematology instruments. Dissertation. Justus-Liebig University, Giessen, Germany). As LYM count in horses showed good agreement, the Mythic 18 can be judged as reliable for counting LYM in horses.

Compared to previous studies of impedance-based hematology instruments, the Mythic 18 showed better agreement with the reference methods for PLT counts for cats, dogs, and horses (Becker M: 2007, A comparative study of seven in-house and two laboratory hematology instruments).⁵ The slight overestimation of canine PLT compared with results of the Sysmex XT-2000iV could be explained by the fact that the Mythic 18 has no sheath flow device. The good results for feline PLT counts were remarkable, particularly as samples with platelet aggregates were included in the calculation. Impedance-based hematology instruments have problems in counting feline PLT accurately.¹⁶ The PLT and RBC sizes in cats often overlap,²⁴ and impedance-based instruments differentiate cells based on their volume. In the present study, more than 53% of the feline samples showed platelet aggregation.




Figure 4. Linearity plot for feline white blood cell count (WBC; **A**), red blood cell count (RBC; **B**), hematocrit value (HCT; **C**), and hemoglobin concentration (HGB; **D**). X-axis: dilution series in %; Y-axis: feline WBC, RBC, HCT, and HGB measured by the Mythic 18 hematology analyzer.^a

Results for the linearity study showed that the Mythic 18 underestimated high WBC values and high HCT values. This is not of severe clinical relevance, as these values were far above the upper reference limit. In the majority of samples analyzed, results of the Mythic 18 would have led to the same clinical interpretation as the results obtained by the reference methods.

Misidentification of platelet aggregates as leucocytes was the most likely cause of lack of detection of leucopenia in 3

feline samples, and identification of erroneous leucocytosis in 10 feline samples. Twenty feline cases with lymphocytopenia were missed with the Mythic 18 due to the presence of platelet aggregates or large platelets. In all cases where the Mythic 18 revealed a lymphocytopenic cell count result, the results were accurate. This does not exclude that occasionally a lymphocytopenia may not be detected if more samples had been tested. For the high degree of misclassification in the LYM count of the dog, no obvious explanation can be provided.




Figure 5. Linearity plot for canine white blood cell count (WBC; **A**), red blood cell count (RBC; **B**), hematocrit value (HCT; **C**), hemoglobin concentration (HGB; **D**), and platelets (PLT; **E**). X-axis: dilution series in %; Y-axis: canine WBC, RBC, HCT, HGB, and PLT measured by the Mythic 18 hematology analyzer.^a




Figure 6. Linearity plot for equine red blood cell count (RBC; **A**), hematocrit value (HCT; **B**), hemoglobin concentration (HGB; **C**), and platelets (PLT; **D**). X-axis: dilution series in %; Y-axis: equine RBC, HCT, HGB, and PLT measured by the Mythic 18 hematology analyzer.^a

The misidentified feline sample with granulopenia showed only a slight difference (7%), and would not have led to a different clinical conclusion. In the 10 cases where the Mythic 18 missed granulocytosis, the cases showed lower total WBC in the Mythic 18 compared to the Sysmex XT-2000iV. In 2 out of the 10 cases, the clinical interpretation would have been different. Platelet aggregates had led to overestimated WBC in 3 feline samples, as described above, which therefore led to false results of

granulocytosis. False-positive granulocytosis and granulocytopenia in equine samples was mainly due to differences in total WBC between the Mythic 18 and the reference instrument.

For RBC, differences in canine samples were all below 10%. For equine HCT values, 4 false-positive samples assuming anemia occurred. The differences were less than 2%, which is attributed to the imprecision of the HCT reading in the capillary tube.

Table 7. Results for carry-over of the Mythic 18 hematology analyzer^a for white blood cell count (WBC), red blood cell count (RBC), hemoglobin concentration (HGB), and platelets (PLT).

Parameter	Species	Sample 1		Sample 2	
		Value	Carry-over (%)	Value	Carry-over (%)
WBC ($\times 10^3/\mu\text{l}$)	Cat	45.8	0.22	9.0	1.11
	Dog	40.1	0.25	5	0
	Horse	9.4	0	5.5	0
RBC ($\times 10^6/\mu\text{l}$)	Cat	3.76	0.26	7.14	0.56
	Dog	5.31	0.19	6.87	0.15
	Horse	7.9	0.25	8.05	0.25
HGB (g/dl)	Cat	4.2	0	10.4	0
	Dog	12.2	0	13.4	0
	Horse	14.4	0	13.8	0
PLT ($\times 10^3/\mu\text{l}$)	Cat	620	0.97	579	0
	Dog	336	0	146	0
	Horse	147	0	147	0

Table 8. Statistically significant changes in blood cell parameters in a cell aging study over 48 hr.

Species	Parameter*	6 hr	10 hr	24 hr	32 hr	48 hr
Cat	MCHC	–	(5.8%) [↓]	(8.5%) [↓]	(10.5%) [↓]	(12.9%) [↓]
	MCV	–	–	(8.6%) [↑]	(9.3%) [↑]	(12.6%) [↑]
	HCT	–	–	–	–	(16.1%) [↑]
	HGB	(3.5%) ^{↑↓}	–	–	–	–
Dog	RBC	–	–	(4.8%) ^{↑↓}	–	–
	MCV	–	–	(5.2%) [↑]	(7.2%) [↑]	(11.2%) [↑]
	MCH	–	–	–	–	(6.4%) [↑]
	MONO%	–	–	–	–	(39.7%) [↑]
Horse	WBC	–	–	(6.6%) [↓]	(18.8%) [↓]	(22.2%) [↓]
	LYM	–	–	(118%) [↑]	(176.6%) [↑]	(184.9) [↑]
	GRAN	–	–	(49.2%) [↓]	(84.4%) [↓]	(90.9%) [↓]
	MCHC	–	–	–	(3.8%) [↓]	(4.2%) [↓]

*MCHC = mean corpuscular hemoglobin concentration; MCV = mean corpuscular volume; HCT = hematocrit value; HGB = hemoglobin concentration; RBC = red blood cell count; MCH = mean corpuscular hemoglobin; MONO% = percentage of monocytes; WBC = white blood cell count; LYM = lymphocytes; GRAN = granulocytes; – = no significant change.

[†]Significant change only one time.

In 19 out of 90 cat samples, the Mythic 18 and the Sysmex XT-2000iV showed thrombocytopenia. Five feline samples with thrombocytopenia were detected only by the Sysmex XT-2000iV. Four out of these 5 feline blood samples demonstrated moderate to severe platelet aggregation in the blood smear. Additionally, in 13 feline samples, the Mythic 18 falsely showed a thrombocytopenia. Platelet aggregation was only found in 1 of these 13 cases, and giant platelets were found in 3 cases. In the remaining cases, no explanation for the detection of false-positive thrombocytopenia can be offered. It has been demonstrated that EDTA anticoagulated blood is prone to build platelet aggregates in cats.¹⁵ High CVs in the equine precision study as well as the narrow range of reference limits may contribute to the high rate of misclassification in equine PLT counts.

In the present study, important abnormalities would have been missed when relying only on the electronic WBC differential of the Mythic 18 (Table 10). Two canine samples showed more than 4 nucleated red blood cells, while 1 feline sample showed 52 nucleated red blood cells. In these samples, WBC results were falsely increased and would have led to different clinical conclusions. One canine sample presented atypical LYM due to an immune-mediated disease. Left shifts and especially degenerative left shifts would have been missed in a remarkable number of blood samples in the canine and feline samples. Foamy cytoplasm of segmented neutrophils has been observed in all 3 species by manual microscopy. This is an important morphological indicator for severe inflammatory disease and toxicity. The presence of reactive LYM is a useful hint to antigenic

Table 9. Clinical relevance of the Mythic 18 hematology analyzer^a with the Sysmex XT-2000iV hematology analyzer^b: results for below, within, and above reference method values.

Parameter*	Species	Agreement			Disagreement		No. of samples
		Below	Within	Above	Below	Above	
WBC	Cat	1/4	64/77	39/48	3	10	129
	Dog	2/2	56/59	60/61	0	3	122
	Horse	9/9	58/61	50/53	0	3	123
LYM (absolute)	Cat	11/31	75/93	3/5	0	18	129
	Dog	6/24	70/90	6/8	5	18	122
	Horse	8/13	87/94	14/16	1	6	123
GRAN (absolute)	Cat	2/3	81/86	30/40	2	3	129
	Dog	0/0	62/64	57/58	0	2	122
	Horse	5/6	106/109	6/8	2	1	123
RBC	Cat	32/35	98/98	7/7	0	0	140
	Dog	45/51	61/63	7/8	0	2	122
	Horse	18/20	81/83	19/20	2	0	123
HCT	Cat	54/58	74/76	3/3	2	0	137
	Dog	59/60	54/58	2/4	4	0	121
	Horse	21/21	86/90	10/12	4	0	123
PLT	Cat	19/24	52/65	1/1	13	0	90
	Dog	13/18	74/82	22/22	0	8	121
	Horse	9/13	69/90	14/19	9	11	122

*WBC = white blood cell count; LYM = lymphocytes; GRAN = granulocytes; RBC = red blood cell count; HCT = hematocrit value; PLT = platelets.

Table 10. Cellular and morphological characteristics associated with misclassification by the Mythic 18 hematology analyzer.^a

Missed abnormalities	Cat	Dog	Horse
Left shift	15/129	29/117	9/123
Normoblasts	12/129	27/122	2/123
Reactive lymphocytes	33/123	10/117	32/123
Cytoplasmatic basophilia and vacuolation	14/123	19/117	21/123
Atypical lymphocytes	0	1	0
Platelet aggregation	48/90	11/121	24/117
Large platelets	7/90	20/121	0

stimulation in the patient.¹⁹ These pathological findings give important information to the clinician and help to improve patient care.

The Mythic 18 was found to perform very well for RBC parameters and total WBC in all investigated species. In cats, it is important to ensure that no platelet aggregates are presented, otherwise WBC and PLT values should be determined by manual methods. The GRAN and LYM counts are accurate in horses. In dogs and cats, absolute GRAN counts are reliable. As with all impedance-based hematological instruments, a microscopic blood smear evaluation is needed to identify platelet aggregates, normoblasts, left shift, cell precursors, and blood parasites and to verify WBC differentiation. Flags for pathological values and reference limits need to be created by the manufacturer of the instrument.

Acknowledgements

The authors thank the technicians of the Clinical Laboratory, Vetsuisse Faculty, University of Zurich, for excellent laboratory assistance, especially Juliette Wälchli. The authors are grateful to Philippe Daire for his excellent support.

Sources and manufacturers

- Orphée SA, Geneva, Switzerland.
- Sysmex Corp., Kobe, Japan.
- Myt-3D, lots B059, B089, B119, Orphée SA, Geneva, Switzerland.
- Mythic 18 Vet M-Pack, Orphée SA, Geneva, Switzerland.
- HemaTek, Siemens, Switzerland.
- Rotina 35 R, Hettrich AG, Switzerland.
- Microsoft Excel® 2007, Microsoft Corp., Redmond, WA.
- Analyse-it Software Ltd., Leeds, United Kingdom.
- Version 3.00 for Windows, GraphPad Software, San Diego, CA.

Declaration of conflicting interests

The authors declared that they had no conflicts of interest with respect to their authorship or the publication of this article.

Funding

This study was funded in part by Orphee SA (Geneva, Switzerland). The authors declared that they received financial support for this project (salary for doctoral candidate AKW). A contract guaranteed independent publication of the results.

References

1. Altman DG, Bland JM: 1983, Measurement in medicine: the analysis of method comparison studies. *Statistician* 32: 307–317.
2. Bablok W, Passing H: 1985, Application of statistical procedures in analytical instrument testing. *J Automat Chem* 7: 74–79.
3. Becker M, Moritz A, Giger U: 2008, Comparative clinical study of canine and feline total blood cell count results with seven in-clinic and two commercial laboratory hematology analyzers. *Vet Clin Pathol* 37:373–384.
4. Boisvert AM, Tvedten HW, Scott MA: 1999, Artifactual effects of hypernatremia and hyponatremia on red cell analytes measured by the Bayer H*1 analyzer. *Vet Clin Pathol* 28:91–96.
5. Deprez P, Bauwens C, Van Schandevijl K, et al.: 2009, Evaluation of the poch-100iV DIFF hematology analyzer for use in horses and cattle. *Vlaams Diergeneeskundig Tijdschrift* 78:105–109.
6. Dewhurst EC, Crawford E, Cue S, et al.: 2003, Analysis of canine and feline haemograms using the VetScan HMT analyzer. *J Small Anim Pract* 44:443–448.
7. International Council for Standardization in Hematology (ICSH): 1994, Guidelines for the evaluation of blood cell analysers including those used for differential leucocyte and reticulocyte counting and cell marker applications. Prepared by the ICSH Expert Panel on Cytometry. *Clin Lab Haematol* 16:157–174.
8. Kjelgaard-Hansen M, Jensen AL: 2006, Is the inherent imprecision of manual leukocyte differential counts acceptable for quantitative purposes? *Vet Clin Pathol* 35:268–270.
9. Knoll JS: 2000, Clinical automated hematology systems. *In: Schalm's veterinary hematology*, ed. Feldman BV, Zinkl JG, Jain NC, 5th ed., pp. 3–11. Lippincott, Williams & Wilkins, Baltimore, MD.
10. Knoll JS, Rowell SL: 1996, Clinical hematology. In-clinic analysis, quality control, reference values, and system selection. *Vet Clin North Am Small Anim Pract* 26:981–1002.
11. Kroll MH, Emancipator K: 1993, A theoretical evaluation of linearity. *Clin Chem* 39:405–413.
12. Laboratory Equipment and Methods Advisory Group: 1969, Recommended scheme for the evaluation of instruments for automatic analysis in the clinical biochemistry laboratory. *J Clin Pathol* 22:278–284.
13. Lilliehöök I, Tvedten H: 2009, Validation of the Sysmex XT-2000iV hematology system for dogs, cats, and horses. I. Erythrocytes, platelets, and total leukocyte counts. *Vet Clin Pathol* 38:163–174.
14. Lilliehöök I, Tvedten H: 2009, Validation of the Sysmex XT-2000iV hematology system for dogs, cats, and horses. II. Differential leukocyte counts. *Vet Clin Pathol* 38:175–182.
15. Moritz A, Hoffmann C: 1997, Thrombozytenzählung bei der Katze [Platelet count in the cat]. *Tierarztl Prax Ausg K Klientiere Heimtiere* 25:695–700. In German.
16. Norman EJ, Barron RC, Nash AS, Clampitt RB: 2001, Prevalence of low automated platelet counts in cats: comparison with prevalence of thrombocytopenia based on blood smear estimation. *Vet Clin Pathol* 30:137–140.
17. Roleff S, Arndt G, Bottema B, et al.: 2007, Clinical evaluation of the CA530-VET hematology analyzer for use in veterinary practice. *Vet Clin Pathol* 36:155–166.
18. Sanzari M, De Toni S, D'Osualdo A, et al.: 1998, Complete analytical and diagnostic performances of the Abbott Cell Dyn 3500. *Panminerva Med* 40:116–125.
19. Stockham SL, Scott MA, eds.: 2008, *Leucocyte. In: Fundamentals of veterinary clinical pathology*, 2nd ed., pp. 53–106. Wiley-Blackwell, Ames, IA.
20. Weiser MG: 1987, Modification and evaluation of a multichannel blood cell counting system for blood analysis in veterinary hematology. *J Am Vet Med Assoc* 190:411–415.
21. Weissenbacher S, Riond B, Hofmann-Lehmann R, Lutz H: 2011, Evaluation of a novel haematology analyser for use with feline blood. *Vet* 187:381–387.
22. Welles EG, Hall AS, Carpenter DM: 2009, Canine complete blood counts: a comparison of four in-office instruments with the ADVIA 120 and manual differential counts. *Vet Clin Pathol* 38:20–29.
23. Wenger-Riggenbach B, Hässig M, Hofmann-Lehmann R, Lutz H: 2006, Evaluation of the LaserCyte: an in-house hematology analyzer for dogs and cats. *Comp Clin Path* 15:117–129.
24. Zelmanovic D, Hetherington EJ: 1998, Automated analysis of feline platelets in whole blood, including platelet count, mean platelet volume, and activation state. *Vet Clin Pathol* 27:2–9.